

Worksheet for tissue processing (up to 12 pieces of tissue)

NAME: _____

1. Tissue processing basket (This should be done the day before processing.)
 1. Take Tissue processing basket(s) out of the acetone jar and air dry.
 2. Rinse the basket(s) a few times with Na cacodylate buffer (0.1M).
2. Epoxy resin (This can be done the day before processing.)
 - Beakers A and B: Mix for ≥ 15 min with the beakers covered with aluminum foil.
 - Beaker A+B: Mix A and B for 15 min, then add DMP-30, followed by mixing for another 15 min. If preparing the day before, wait until the next day to mix DMP-30.

	For Day 1			For Day 2		
	A	B	A+B	A	B	A+B
LX112	7.00 g	18.00 g	-	14.00 g	36.00 g	-
DDSA	-	-	-	-	-	-
NMA	-	-	-	-	-	-
A	-	-	13.5 g	-	-	27 g
B	-	-	31.5 g	-	-	63 g
DMP-30	-	-	0.63 g	-	-	1.26 g

NOTE: The amounts of DDSA and NMA must be adjusted according to the WPE value indicated on the LX-112 bottle. Refer to the table in the full protocol and record their amounts in this table. Also record WPE of LX-112 here: _____

3. Reagents for osmium fixation, dehydration, and *en bloc* UA staining

- In 50-ml conical tubes, add:
 - ddH₂O: 25 ml for setup + 20 ml for processing
 - Na cacodylate buffer (SCB; 0.2M): 12 ml for setup
 - Na cacodylate buffer (0.1M): ~150 ml total for processing
- In the labeled scintillation vials, prepare the following (Mix well by shaking after adding each reagent):

Reagents (add in this order ↓)	Reduced Osmium*			1% OsO ₄ **	EtOH 50%	EtOH + UA				EtOH 100%
	1	2	KFeCN			50%	70%	90%	100%	
ddH ₂ O	-	-	5 ml	3 ml	5 ml	5 ml	3 ml	1 ml	-	-
KFeCN	-	-	0.3 g	-	-	-	-	-	-	-
0.2M SCB	3 ml	3 ml	-	6 ml	-	-	-	-	-	-
4% OsO ₄	3 ml	-	-	3 ml	-	-	-	-	-	-
EtOH	-	-	-	-	5 ml	-	2 ml	4 ml	5 ml	10 ml
2% UA	-	-	-	-	-	5 ml	5 ml	5 ml	5 ml	-

*Reduced osmium:

1. Add 3 ml of KFeCN solution to vial 2.
2. Place vials 1 and 2 on ice.
3. When ready, pour vials 1 and 2 back and forth to mix well before adding to tissue.

**1% OsO₄: Keep vial on ice until use.

4. Reagents for infiltration

- Prepare ~2 ml per tissue of the following:
 - 1:1 = EtOH : propylene oxide (PO)
 - 1:2 = EtOH : PO
 - 100% PO
 - 1:1 = PO : Resin
 - 1:2 = PO : Resin

5. Microwave

1. Fill load cooler with fresh RO water and eliminate air bubbles. Make sure pump is on.
2. Set power level to 1 (175W).
3. Run microwave for 2 min to warm up magnetron.


Worksheet for tissue processing (up to 12 pieces of tissue)

NAME: _____

6. Processing (4 days)

DAY 1 (Date: _____)

- Place a piece of absorbent paper on work surface in the fume hood. Prepare new solid waste bags.
- Embed vibrasliced tissue in 7-9% agarose, and place them into processing basket filled with Na cacodylate buffer (0.1M).



Description of tissue (e.g., animal ID, vibraslice ID, experiment description, genotype, etc.):

Worksheet for tissue processing (up to 12 pieces of tissue)

NAME: _____

Waste: Aldehydes-Caco Osmium-KFeCN

DAY 1 continued (Date: _____)

- Na cacodylate buffer (0.1M) rinses, 5 min × 5 _____
- Reduced osmium, 5 min _____
- Na cacodylate buffer (0.1M), 5 min × 5 _____
- Microwave power setting 1 (175W), under vacuum
 - 1% OsO₄, 1 min ON → 1min OFF → 1 min ON _____
 - Cool the tissue to about 15°C.
 - 1% OsO₄, 1 min ON → 1min OFF → 1 min ON _____
- Na cacodylate buffer (0.1M) rinses, 2 min × 5 _____
- ddH₂O rinses, 2 min × 2 _____

UA-EtOH


- 50% EtOH, 5 min _____
- Microwave power setting 2 (250W), no vacuum
 - 50% EtOH + UA, 40 s _____
 - 70% EtOH + UA, 40 s _____
 - 90% EtOH + UA, 40 s _____
 - 100% EtOH + UA, 40 s _____
 - 100% EtOH, 40 s _____ [Turn off the pump on microwave.]

Flammable PO-Resin

- Transfer tissue into shell vials filled with 1:1 = EtOH : PO.
- Infiltration (on rotator)
 - 1:1 = EtOH : PO, 10 min _____
 - 1:2 = EtOH : PO, 10 min _____
 - 100% PO, 15 min × 2 _____
 - 1:1 = PO : Resin, 1 hr _____
 - 1:2 = PO : Resin, overnight _____

DAY 2 (Date: _____)

- 100% Resin, 1 hr × 3 _____
- Transfer tissue into mold, fill with resin, and place in 60°C oven for 48 hrs. Current time is: _____:_____.



Worksheet for tissue processing (up to 12 pieces of tissue)

NAME: _____

DAY 2 continued (Date: _____)

- Clean up after yourself.
 - Excess liquid resin: Collect into a disposable beaker along with other resin-contaminated supplies (all vials must be uncapped), then place in 60°C oven for 48 hrs.
 - Solid waste: All vials must be uncapped. Place all UA-contaminated waste (incl. absorbent paper) in a dedicated bag. Other solid waste should be placed in a “No UA” waste bag.
 - 12- or 24-well plates with vibraslices: Discard all liquid into “Aldehydes-Caco” waste bottle. Then place the plates into “No UA” solid waste bag.

DAY 3 (Date: _____)

- Continue resin polymerization. You may check and adjust the positioning of tissue at this point.

DAY 4 (Date: _____ , Time: _____)

- Take the resin-embedded tissue blocks out of oven. Remove the blocks out of the mold immediately and store in the cardboard pill boxes labeled with block identification numbers.
- Take solidified resin waste into “Solid Resin” waste drum.
- Examine the blocks under the stereomicroscope and record their images before trimming and cutting.
- Scan this note into a pdf file.