


Performing a Retro-reconstruction on the Siemens Vida (SyngoMR XA)

To reconstruct DICOM images from raw data

1) Open the **Job View** tool (Figure 1), by clicking on the **Briefcase** icon in the upper right-hand corner of the **Data Management** window. You can also click the computer icon in the lower left-hand corner of the window.



Figure 1: Job View Tool

- 2) In the **Job View** tool, select the **MRI Reconstructions** tab: The tool will display raw-data (Figure 2) that remains in the system cache. These are listed by measurement number (**MeasUID**), **Patient name** and **Scan Protocol** name.
- 3) Select the measurement / scan you wish to reconstruct.
- 4) Under the **Action** column click the curly arrow button  to start the reconstruction of the scan data; any dependent measurements (e.g coil calibrations) will also be reconstructed.




Figure 2: MRI Reconstructions tab

nb/ depending on the scan-type, a reconstruction can take several minutes to perform, particularly for SMS-EPI sequences with a large # of repetitions.