


Configure VLC to output to MRI Speakers

This document describes how to change settings in VLC media player to output audio to the MRI scanner speakers rather than headphones. Typically this is used to play MRI noise to the speakers.

Open the file you wish to play.

MRI recording is saved at C:
\\Users\\Public\\Music\\sfx\\MRI 3T Scanner
Noise

In VLC, click Audio > Audio Device >
Line Out (Scarlett 6i6 USB)



Adjust the volume to the MRI using the
large silver knob on the red Focusrite
Box on the desk.



Reviewed 10/2/17 - Cox